

1. Collect feces specimen in the cardboard container provided. Do not contaminate with urine or toilet tissue.
2. Label the 4-ounce plastic container (not the lid) with the patient's name, date of birth, and the date and time of the collection.
3. Transfer a sizable amount of feces, no smaller than the size of a "pea," to the 4-ounce container, tightly secure the lid and place in the biohazard bag. Freeze specimen immediately and transport to the lab frozen.
4. The specimen must be delivered to the lab within 24 hours of collection.

1. Recoja la muestra de materia fecal en la recipiente que se proveé. No lo contamine con orina o papel higienico.
2. Etiquetar el envase de plástico de 4 onzas (no la tapa) con el nombre del paciente, la fecha de nacimiento y la fecha y hora de la recolección.
3. Transferir una cantidad considerable de la materia fecal no mas pequeño que el tamaño de un guisante dentro del contenedor de plástico de cuatro onzas. Asegure firmemente la tapa y colóquela en la bolsa bio-peligro. Congelar el espécimen de inmediato y transportarlo al laboratorio.
4. El espécimen debe ser entregado al laboratorio dentro de veinticuatro horas de la colección.