



LAB ALERT

THE MEDICAL FOUNDATION

Blood Culture Changes

The Medical Foundation is converting to a new blood culture system.

Bottles and collection devices are changing.

**New culture bottles and their compatible transfer devices will be available
January 1, 2017.**

	TMF Stock No.	Unit of Issue
BacT/Alert® FA Plus Aerobic / Fungus	9420	Each
BacT/Alert® FN Plus Anaerobic	9423	Each
BacT/Alert® PF Plus Aerobic Pediatric	9422	Each

Jelco® Saf-T Holder® and safety insert – sterile Smiths Medical Compatible with female Luer collection devices	3562	Each
--	------	------

Monoject™ Transfer Device and safety insert –sterile Covidien™ Monoject™ Compatible with male Luer collection devices	3563	Each
---	------	------

Revised supply order forms are available from Client Services

Conversion to the new bottles must be completed by March 31, 2017.

Unused BD bottles may be returned if not depleted prior to conversion.

Instructive videos demonstrating the new devices may be viewed online at:

[BacT-Alert-direct-draw.html](https://www.biomerieuxuniversity.com/bact-alert-direct-draw.html)

<https://www.biomerieuxuniversity.com/bact-alert-direct-draw.html>

[BacT-Alert-syringe-draw.html](https://www.biomerieuxuniversity.com/bact-alert-syringe-draw.html)

<https://www.biomerieuxuniversity.com/bact-alert-syringe-draw.html>

THE MEDICAL FOUNDATION

530 North Lafayette Boulevard • South Bend, Indiana 46601
(574) 234-4176 • Elkhart (574) 293-8441 • (800) 544-0925
Joyce L. Simpson, M.D. • *Medical Director*

The Medical Foundation is converting to a new blood culture system

Culture supplies are changing
 Conversion to the new bottles must be completed by
March 31, 2017



AEROBIC / Fungus
 NEW Order No.
 9420



ANAEROBIC
 NEW Order No.
 9423



PEDIATRIC
 NEW Order No.
 9422

These bottles will require new collection devices:



Male Luer Hub with safety insert
 Compatible with female Luer
 collection devices

NEW Order No.
 3562



Female Luer Hub with safety insert
 Compatible with male Luer
 collection devices

NEW Order No.
 3563

Revised supply order forms are available from Client Services
 Unused BD bottles may be returned if not depleted prior to conversion

THE MEDICAL FOUNDATION

530 North Lafayette Boulevard • South Bend, Indiana 46601
 (574) 234-4176 • Elkhart (574) 293-8441 • (800) 544-0925
 Joyce Simpson, M.D. • Medical Director