

**The tests listed in this table will be DELETED from our catalog.  
Replacement or alternate tests are noted if applicable.**

<b>Order Code</b>	<b>Effective Date</b>	<b>Deleted Test</b>	<b>Replacement / Alternate Test</b>
35342	2.21.2017	Enterovirus Antibody Profile (Care Set)	Enterovirus Antibodies Panel (47564) – See January 2017 notice.
43656	Immediate	Sin Nombre Virus (SNV) IgG Confirmation, IBL	N/A
43657	Immediate	Sin Nombre Virus (SNV) IgM Confirmation, ELISA	N/A

**NEW TEST 28437 Allergens, Xolair Perennial Allergy Panel, IgE**

Available Immediately

Performed by The Medical Foundation

Specimen information sheet enclosed.

**36350 Chlamydia trachomatis and Neisseria gonorrhoeae by NAA, Swab**

Effective Immediately

**Rectal** and **oral pharyngeal** swabs are now acceptable for this test based on in-house validation by The Medical Foundation. Collect using the **cobas®** PCR Female Swab Sample Kit.

Note, these sources are off-label.

This change also applies to:

- 36351 Chlamydia trachomatis by NAA, Swab
- 36352 Neisseria gonorrhoeae by NAA, Swab

**23057 Electrophoresis, Protein, with Interpretation, Serum**

Effective Immediately

Test reactivated

Test Name changed. Previously listed as “Electrophoresis, Protein, Serum”

Serum protein fractionation by capillary electrophoresis is performed on all cases.

This panel does not include a Monoclonal Protein ID (Immunofixation or Immunotyping) assay.

Reports include pathologist interpretation. An additional professional fee applies. Please provide a face sheet with the patient’s billing information.

Specimen information sheet enclosed.

**47094 FLT3 Mutation Detection, PCR**

Effective Immediately

Changes to specimen type, minimum volume, stability, and storage/transport.

Specimen information sheet enclosed.

**44092 Fluoxetine and Metabolite, Serum or Plasma**

Effective Immediately

Performing Lab changed to NMS Labs

Specimen information sheet enclosed.

**43614 Hantavirus Antibodies, IgG and IgM by ELISA**

Effective Immediately

Reflex to Sin Nombre Virus IgG and/or IgM confirmation will no longer be performed.

---

**NEW TEST    28011    Hepatitis C Virus Antibody with Reflex to HCV RNA,  
Quantitative RT-PCR and Genotyping**

Available Immediately

Performed by The Medical Foundation

Specimen information sheet enclosed.

If HCV Antibody result is REACTIVE, then HCV RNA, Quantitative RT-PCR will be added.

If HCV viral load is 1000 IU/mL or greater, then HCV Genotyping will be added.

The following reflex tests may apply:

- 36021 – Hepatitis C Virus RNA, Quantitative RT-PCR with Reflex to Genotyping
- 36123 – Hepatitis C Virus Genotyping by RT-PCR

---

**36004            Herpes Simplex Virus Type 1 and 2, Qualitative PCR**

Effective Immediately

Blood and Pap vial specimen types have been removed from the collection information for this test. Use this test for swab specimens of genital tract, skin, or mucous membranes. Refer to specific tests for other specimen types:

- Blood – Test code 36012
- CSF – Test code 36013
- Pap vial – Test code 36014

Specimen information sheet enclosed.

---

**NEW TEST    36012    Herpes Simplex Virus Type 1 and 2, Qualitative PCR, Blood**

Available Immediately

Performed by The Medical Foundation

Specimen information sheet enclosed.

---

**45708            Ibuprofen, Serum or Plasma**

Effective Immediately

Performing Lab changed to NMS Labs

Specimen information sheet enclosed.

---